

BLOOD GLUCOSE MONITORING SYSTEM



±15mg/dL & ±15%



LCD Backlight



Strip Ejection



Strip Guiding Light

5 Sec, 0.5 μ L

Display Screen

Strip-Ejection



FC-G168

FEATURES

- Strip port indicator
- Large and backlit display
- Before/After meal (AC/PC) presetable
- Connecting the FaCare app on smartphones or tablets
- Support data backup to Cloud Server

SPECIFICATIONS

Enzyme Type	GDH-FAD
Sample Size	0.5 μ L
Reaction Time	5 seconds
Measurement Range	20 ~ 600 mg/dL (1.1 ~ 33.3 mmol/L)
Hematocrit Range	35% - 60%
Precision	SD < 5 mg/dL (0.278 mmol/L) at < 100 mg/dL (5.55 mmol/L); CV < 5% at 100 mg/dL (5.55 mmol/L)
Accuracy	±15 mg/dL if < 100 mg/dL; ±15% if ≥ 100 mg/dL
Ketone Warning	Yes
Communication	Bluetooth
Power Source	AAA x 2
Memory Capacity	1000 sets
Day Average	7, 14, 21, 28, 60, 90 days
Daily Alarm	4 daily alarms
Dimension	96 (L) x 61 (W) x 26 (H)mm
Weight	67.2 g (without battery)
Operating Condition	10°C (50°F) ~ 40°C (104°F), below 85% R.H.
Storage Condition	-20° C (-4°F) ~ +60°C (140°F) (Meter) 2°C (35.6°F) ~ 32°C (89.6°F) (Strip)



Q FaCare



FaCare

Android

IOS

FaCare⁺

FACARE INTERNATIONAL MEDICAL TECHNOLOGY JSC.

Ha Noi office: 9th Floor, N01T3 Building, Doan Ngoai Giao Zone, Xuan Tao Ward, Bac Tu Liem District, Ha Noi
Tel: (+84-24) 6295 5592 - Fax: (+84-24) 37914142
Ho Chi Minh office: B24 Building, Nguyen Anh Thu Project, Quarter 3, Tan Chanh Hiep Ward, 12 District, Ho Chi Minh City
Web: www.facare.com.vn - Email: info@facare.com.vn

Family doctor at your home!

Hotline: 0932 372 099 - 0906 209 266